

Low Vision Aid

**Prismatic Spectacles**

Model:	Prismatic Spectacles
PS68:	+6D spherical, 8D prism
PS810:	+8D spherical, 10 D prism
PS1012:	+10 D spherical, 12 D prism
PS1214:	+12 D spherical, 14 D prism

All prisms are base inwards.

Protective and anti-reflective BBAR coatings are an optional choice.

Contact Us



TARUN ENTERPRISES

[ADDRESS](#)

8/8 Strachy Road, Prayagraj, Uttar Pradesh, INDIA – 211001

Website: www.optitecheyecare.com | Email: info@optitecheyecare.com | Phone: +91-9335154556